

## CRYOSAVE I

DMEM/F12, 5% DMSO, Human Serum Albumin (USP), Amino Acids, Electrolytes

**Catalog number:** 136 (100 mL), 137 (500 mL), 184 (10 mL), 185 (50 mL)

**Size:** 10 mL, 50 mL, 100 mL and 500 mL

**Intended Use:**

For research purposes only.

Cryosave I is a cryo-preservation medium for freezing cells, ranging from mesenchymal stem cells, hepatocytes, other cell types to tissue sample, at -86°C, and -196°C.

**Summary and Explanation:**

Cryosave I is a ready-to-use and complete medium with no further additives required. The formula was improved with a reduction of DMSO (10%) and xeno-free. In Cryosave I, the fetal bovine serum is substituted for human serum albumin, minimizing the risk of xeno-immunization and zoonotic transmission.

Cryosave I maintains high cell viability on post-thawing, which is ideal for scientific research and research on technology transfer.

**Known Applications:**

Cryosave I demonstrates a cryopreservation efficiency, which results in a high survival rate of thawed cells, more than 85%, for mesenchymal stem cells derived from various sources, including adipose, umbilical cord tissue, and bone marrow.

**Reagents Provided:**

- 10 mL Cryosave I or
- 50 mL Cryosave I or
- 100 mL Cryosave I or
- 500 mL Cryosave I

**Reconstitution, Dilution, and Mixing:**

Directly resuspend cells in Cryosave I.

Cryosave I is provided at an 1x concentration to be added directly to cells without dilution. Dilution or mixing is not necessary.

**Materials and Reagents Required But Not Provided:**

Not applicable

**Storage and Stability:**

Stored at -20 to 4 °C.

Shelf life at 12 months.

**Instructions for Use:**

*Cell Freezing:*

1. After harvesting, centrifuge cells at 1.500 rpm for 5 minutes.
2. Remove supernatant and resuspend in Washing Buffer for cell washing.
3. Re-centrifuge and collect cell pellets.
4. Resuspend cell pellet in cold Cryosave I (2-8 °C) at a cell density of 1-2 .10<sup>6</sup> cells per mL. Mix thoroughly to achieve a homogeneous cell suspension.
5. Aliquot cell suspension into cryogenic vials.

*Cooling:*

- **6.1.** Freeze vials at -20°C in 120 minutes; then transfer vials to -80°C or -196°C (liquid nitrogen for long-term storage). *or*
- **6.2.** Put cryogenic vials into cooling box for -1 °C/min (for example, Mr. Frosty), transfer boxes to -86°C overnight; finally transfer the vials to -196°C (liquid nitrogen for long-term storage).

*Thawing Cells:*

1. Remove the cryogenic vials from storage.
2. Place vials into a 37°C water bath for 1-2 min.
3. Then transfer the vial into a laminar hood.
4. Transfer thawed cells into centrifuge tube containing Thawbest or media. Appropriate amount of Thawbest (ratio 1:4 Cryosave : Thawbest) or pre-warmed growth media into centrifuge tube (ratio 1:4 Cryosave : fresh cultured medium).
5. Centrifuge the cell suspension at 1.500 rpm for 5 minutes.

6. Discard the supernatant without disturbing the cell palletete.
7. Resuspend cells in Washing Buffer or growth medium according to the experimental design.

**Limitations:**

Non-injectable and non-transfusable.

**Quality Control:**

- pH: 7.0 - 8.0
- Bacterial and Fungal Contamination: Negative
- Mycoplasma: Negative.
- Endotoxin: < 1.0 EU/mL
- Colour: Yellow-Red
- Cell viability after thawing: ≥ 85%

**Precautions:**

Rapid or uncontrolled freezing may induce intracellular crystallization, which could affect the viability of thawed cells.

Post thawing, Cryosave I should be discarded, and cells should be washed before proceeding.









Personal protective equipment is required.

**Troubleshooting:**

Not applicable

**Explanation of symbols and warnings**

The symbols on produce labels are explained below:

			
Use By:	Batch code	Keep away from light	Catalog number
			
Temperature Limitation	Consult instructions for use	Caution, consult accompanying documents	Sterilized using aseptic processing techniques

**Related products**

Products	Catalog No.
<b>Washing Buffer</b>	
100 mL	149
500 mL	150
4 L (Bag)	151
<b>Deattachment</b>	
100 mL	120
500 mL	121
200 mL (Bag)	122
<b>ThawBest</b>	
100 mL	142
500 mL	143

To purchase other products, please visit:

<http://biomedmart.org>

For further information: Please contact us at:

[contact@sci.edu.vn](mailto:contact@sci.edu.vn)

[sales@sci.edu.vn](mailto:sales@sci.edu.vn)